



**CERTIFICATE OF ANALYSIS  
Anti-TPO ELISA (96 TESTS)**

**PRODUCT CODE** BXE0896A  
**PRODUCT** ANTI-TPO ELISA  
**BATCH NUMBER** ATPO-1706-1  
**DATE OF MANUFACTURE:** 2017 / 06  
**EXPIRY DATE** 2018 / 05

**RESULTS**

STANDARDS	ABS.	QC TEST	RESULTS	SPECIFICATIONS	OUTCOME
"A" 0.0 IU/ML	0.027	ABS(A)	0.027	<0.08	PASS
"B" 25 IU/ML	0.313	ABS(B/A)	11.59	>2.00	PASS
"C" 50 IU/ML	0.523	ABS(C/B)	1.67	>1.25	PASS
"D" 100 IU/ML	0.925	ABS(D/C)	1.77	>1.50	PASS
"E" 250 IU/ML	1.790	ABS(E/D)	1.94	>1.50	PASS
"F" 500 IU/ML	2.529	ABS(F/E)	1.41	>1.20	PASS
100 % OF "F"	2.529	ABS(F)	2.529	>1.80	PASS

This lot of ANTI-TPO ELISA kit has been tested by quality control department and meets all final Q.C specifications. This product is released for packing and is suitable for use until the stated expiry dates.

PASS / FAIL: **Pass**

SIGNATURE: *Cherith Morris* (QC Manager)

DATE: **2017/06**