

CERTIFICATE OF ANALYSIS

FSH EIA TEST KIT

Catalog No. BC-1029		Lot No. RN-55648		Exp. Date: 01-31-2018	
QC Reference No. 25337			Manufactured Date: 07-26-2016		
Component		Part No.	Lot No.	Exp. Date	
1. FSH Ab Coated Wells	1 ea	MB 20015	RN-55016	04-2019	
2. Standards	1 set	MB 30015	RN-55141	05-2019	
3. Conjugate Reagent	13 ml	MB 40015	RN-55573	01-2018	
4. TMB Reagent	11 ml	MB 10051	RN-55525	07-2018	
5. Stop Solution	11 ml	MB 10074	RN-55293	05-2019	
7. Bio-Rad Controls (not included)		Bio-Rad P/N: 370	40320	05-2018	

STANDARD CURVE

Concentration (mIU/ml)	O.D. 450 nm
0	0.058
5	0.104
15	0.234
50	0.888
100	1.636
200	3.097

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

CONTROL VALUES

Control No.	Mean (mIU/ml)	Range
1	10.4	7.22 - 13.0
2	23.1	16.6 - 31.6
3	53.8	37.0 - 73.9

Pass Reject
Signature  Date 08-08-2016
Quality Control