

CERTIFICATE OF ANALYSIS

FSH EIA TEST KIT

Catalog No. BC-1029		Lot No. RN-52473		Exp. Date: 07-2016	
QC Reference No. 24267			Manufactured Date: 01-16-2015		
Component		Part No.	Lot No.	Exp. Date	
1. FSH Ab Coated Wells	1 ea	MB 20015	RN-52249	11-2017	
2. Standards	1 set	MB 30015	RN-52452	01-2018	
3. Conjugate Reagent	13 ml	MB 40015	RN-52453	07-2016	
4. TMB Reagent	11 ml	MB 10051	RN-52387	01-2017	
5. Stop Solution	11 ml	MB 10074	RN-52333	12-2017	
7. Bio-Rad Controls (not included)		Bio-Rad P/N: 370	40290	04-2016	

STANDARD CURVE

Concentration (mIU/ml)	O.D. 450 nm
0	0.050
5	0.111
15	0.251
50	0.752
100	1.554
200	2.833

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

CONTROL VALUES

Control No.	Mean (mIU/ml)	Range
1	6.39	4.37 - 9.41
2	23.3	20.0 - 30.8
3	55.3	49.3 - 76.3

Pass Reject

Signature James Ho Date 01-16-2015
Quality Control