

CERTIFICATE OF ANALYSIS

**FSH EIA TEST KIT**

<b>Catalog No. BC-1029</b>		<b>Lot No. RN-54772</b>		<b>Exp. Date: 08-31-2017</b>	
<b>QC Reference No. 25070</b>			<b>Manufactured Date: 02-26-2016</b>		
<b>Component</b>		<b>Part No.</b>	<b>Lot No.</b>	<b>Exp. Date</b>	
1. FSH Ab Coated Wells	1 ea	MB 20015	RN-54304	11-2018	
2. Standards	1 set	MB 30015	RN-54743	02-2019	
3. Conjugate Reagent	13 ml	MB 40015	RN-54627	08-2017	
4. TMB Reagent	11 ml	MB 10051	RN-54659	02-2018	
5. Stop Solution	11 ml	MB 10074	RN-54661	02-2019	
7. Bio-Rad Controls (not included)		Bio-Rad P/N: 370	40290	04-2016	

**STANDARD CURVE**

<b>Concentration (mIU/ml)</b>	<b>O.D. 450 nm</b>
0	0.058
5	0.121
15	0.254
50	0.850
100	1.631
200	2.808

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

**CONTROL VALUES**

<b>Control No.</b>	<b>Mean (mIU/ml)</b>	<b>Range</b>
1	6.14	4.37 - 9.41
2	22.3	20.0 - 30.8
3	56.7	49.3 - 76.3

Pass  Reject

Signature Jenna N. Date 02-29-2016  
Quality Control