

CERTIFICATE OF ANALYSIS

Free T4 EIA TEST KIT

Catalog No. BC-1008	Lot No. RN-55592	Exp. Date: 10-31-2017		
QC Reference No. 25335	Manufactured Date: 07-25-2016			
Component		Part No.	Lot No.	Exp. Date
1. fT4 Ab Coated Wells	1 ea	MB 20051	E-41115	11-2018
2. Standards, Lyq.	1 set	MB 30051	12(01-06)136	03-2019
3. Enzyme Conjugate Reagent	10.5 ml	MB 40051	1210136	03-2019
4. Color Reagent A	13 ml	MB 10024	RN-55065	10-2017
5. Color Reagent B	13 ml	MB 10028	RN-55066	10-2017
6. Stop Solution	10 ml	MB 10062	RN-55041	04-2019
7. Bio-Rad Controls (not included)		Bio-Rad P/N: 370	40320	05-2018

STANDARD CURVE

Concentration (µIU/ml)	O.D. 450 nm
0.00	3.053
0.38	2.724
1	1.953
2	1.178
3.5	0.671
7.5	0.325

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

CONTROL VALUES

Control No.	Mean (µIU/ml)	Range
1	0.741	0.263 - 1.27
2	1.96	1.38 - 2.67
3	3.44	2.40 - 4.65

Pass Reject

Signature *[Signature]* Date 07-25-2016

Quality Control