

CERTIFICATE OF ANALYSIS

Free T4 EIA TEST KIT

Catalog No. BC-1008	Lot No. RN-56142	Exp. Date: 02-28-2018		
QC Reference No. 25533		Manufactured Date: 11-04-2016		
Component		Part No.	Lot No.	Exp. Date
1. fT4 Ab Coated Wells	1 ea	MB 20051	E-4196	09-2019
2. Standards, Lyq.	1 set	MB 30051	12(01-06)136	03-2019
3. Enzyme Conjugate Reagent	10.5 ml	MB 40051	1210176	07-2019
4. Color Reagent A	13 ml	MB 10024	RN-55643	02-2018
5. Color Reagent B	13 ml	MB 10028	RN-55644	02-2018
6. Stop Solution	10 ml	MB 10062	RN-55041	04-2019
7. Bio-Rad Controls (not included)		Bio-Rad P/N: 370	40320	05-2018


STANDARD CURVE

Concentration (μIU/ml)	O.D. 450 nm
0.00	3.087
0.38	2.806
1	2.226
2	1.415
3.5	0.847
7.5	0.348

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

CONTROL VALUES

Control No.	Mean (μIU/ml)	Range
1	0.678	0.263 - 1.27
2	1.88	1.38 - 2.67
3	3.27	1.83 - 3.65

Pass Reject
Signature  Date 11-07-2016
Quality Control