

CERTIFICATE OF ANALYSIS

Free T4 EIA TEST KIT

Catalog No. BC-1008	Lot No. RN-54725	Exp. Date: 04-30-2017		
QC Reference No. 25052		Manufactured Date: 02-17-2016		
Component		Part No.	Lot No.	Exp. Date
1. fT4 Ab Coated Wells	1 ea	MB 20051	E-41115	11-2018
2. Standards, Lyq.	1 set	MB 30051	12(01-06)195	09-2018
3. Enzyme Conjugate Reagent	10.5 ml	MB 40051	1210185	08-2018
4. Color Reagent A	13 ml	MB 10024	RN-54189	04-2017
5. Color Reagent B	13 ml	MB 10028	RN-54190	04-2017
6. Stop Solution	10 ml	MB 10062	RN-54306	11-2018
7. Bio-Rad Controls (not included)		Bio-Rad P/N: 370	40290	04-2016

STANDARD CURVE

Concentration (μIU/ml)	O.D. 450 nm
0.00	2.637
0.30	2.396
1.05	1.870
2.15	1.121
3.60	0.618
8.05	0.203

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

CONTROL VALUES

Control No.	Mean (μIU/ml)	Range
1	0.866	0.550 - 1.11
2	2.27	1.28 - 2.51
3	3.12	1.83 - 3.65

Pass Reject

Signature James Rho Date 03-02-2016

Quality Control