

CERTIFICATE OF ANALYSIS

LH EIA TEST KIT

Catalog No. BC-1031	Lot No. RN-55977	Exp. Date: 04-30-2018		
QC Reference No. 25458		Manufactured Date: 09-23-2016		
Component		Part No.	Lot No.	Exp. Date
1. LH Ab Coated Wells	1 ea	MB 20016	RN-55874	09-2019
2. Standards	1 set	MB 30016	RN-55142	05-2019
3. Conjugate Reagent	13 ml	MB 40016	RN-55040E	04-2018
4. TMB Reagent	11 ml	MB 10051	RN-55784	08-2018
5. Stop Solution	11 ml	MB 10074	RN-55738	08-2019
6. Bio-Rad Controls (not included)		Bio-Rad P/N: 370	40320	05-2018

STANDARD CURVE

Concentration (mIU/ml)	O.D. 450 nm
0	0.048
5	0.136
15	0.312
50	0.981
100	1.721
200	2.761

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

CONTROL VALUES

Control No.	Mean (mIU/ml)	Range
1	2.48	1.08 - 3.89
2	19.7	16.0 - 26.3
3	51.1	40.6 - 67.6

Pass Reject

Signature Jenna N Date 11-04-2016
Quality Control