

CERTIFICATE OF ANALYSIS

LH EIA TEST KIT

Catalog No. BC-1031	Lot No. RN-52883	Exp. Date: 08-2016		
QC Reference No. 24406		Manufactured Date: 03-12-2015		
Component		Part No.	Lot No.	Exp. Date
1. LH Ab Coated Wells	1 ea	MB 20016	RN-52755	02-2018
2. Standards	1 set	MB 30016	RN-52599	03-2018
3. Conjugate Reagent	13 ml	MB 40016	RN-52618	08-2016
4. TMB Reagent	11 ml	MB 10051	RN-52723	02-2017
5. Stop Solution	11 ml	MB 10074	RN-52797	02-2018
6. Bio-Rad Controls (not included)		Bio-Rad P/N: 370	40290	04-2016

STANDARD CURVE

Concentration (mIU/ml)	O.D. 450 nm
0	0.047
5	0.144
15	0.332
50	0.998
100	1.732
200	2.819

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

CONTROL VALUES

Control No.	Mean (mIU/ml)	Range
1	1.96	1.00 - 3.18
2	23.4	17.4 - 28.1
3	59.3	40.2 - 73.2

Pass Reject

Signature James M Date 03-30-2015
Quality Control