

CERTIFICATE OF ANALYSIS

**Testosterone EIA TEST KIT**

<b>Catalog No. BC-1115</b>	<b>Lot No. RN-53163</b>	<b>Exp. Date: 11-2016</b>		
<b>QC Reference No. 24373</b>		<b>Manufactured Date: 05-06-2015</b>		
<b>Component</b>		<b>Part No.</b>	<b>Lot No.</b>	<b>Exp. Date</b>
1. Goat Anti-Rabbit Ab Coated Wells	1 ea	MB 20049	RN-52879	03-2018
2. Standards Set	1 set	MB 30049	RN-50731	04-2017
3. Controls (1&2)	1 set	MB 30049C	RN-50732	04-2017
4. Conjugate Reagent	12 ml	MB 40049	RN-53142	11-2016
5. Antibody Reagent	7 ml	MB 70049	RN-53143	05-2017
6. TMB Reagent	11 ml	MB 10051	RN-53150	05-2017
7. Stop Solution (1NHCl)	11 ml	MB 10074	RN-53131	04-2018


**STANDARD CURVE**

<b>Concentration (ng/ml)</b>	<b>O.D. 450 nm</b>
0	2.572
0.1	2.246
0.5	1.641
2.0	1.235
6.0	1.020
18.0	0.638

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

**CONTROL VALUES**

<b>Control No</b>	<b>Mean (ng/ml)</b>	<b>Range (pg/ml)</b>
1	0.466	0.326 - 0.725
2	9.73	7.89 - 12.6

Pass  Reject   
Signature  Date 05-29-2015  
Quality Control