

CERTIFICATE OF ANALYSIS

TSH EIA TEST KIT

Catalog No. BC-1001	Lot No. RN-44382	Exp. Date: 05-2013		
QC Reference No. 19227		Manufactured Date: 12-16-2011		
Component		Part No.	Lot No.	Exp. Date
1. TSH Ab Coated Wells	1 ea	MB 20001	RN-44361	12-2014
2. Standards, Lyo.	1 set	MB 30001	RN-44338	12-2014
3. Conjugate Reagent	13 ml	MB 40001	RN-44219	05-2013
4. TMB Reagent	11 ml	MB10051	RN-44332	12-2013
5. Stop Solution	11 ml	MB10074	RN-44303	12-2014
6. Bio-Rad Controls (not included)		Bio-Rad P/N: 370	40230	05-2012

STANDARD CURVE

Concentration (μIU/ml)	O.D. 450 nm
0	0.052
0.5	0.154
2.0	0.424
5.0	0.895
10.0	1.511
25.0	2.763

This Standard Curve was obtained on the date of performance testing of the above components. It is for the purpose of illustration only. Each laboratory should establish its own standard curve.

CONTROL VALUES

Control No.	Mean (μIU/ml)	Range
1	0.388	0.207 - 0.961
2	5.86	3.99 - 7.75
3	>25.0	22.4 - >25.0

Pass Reject

Signature James H. Date 12-19-2011
Quality Control