



CERTIFICATE OF ANALYSIS
FINISHED KIT QUALITY CONTROL

PRODUCT NAME: **RAPID PLASMA REAGIN (RPR)** MANUFACTURE DATE: 05/04/2011
 CATALOG#: RPR-500
 LOT#: 32355 EXP. DATE: 04/2013

COMPONENTS: 1. RPR ANTIGEN LOT #: 730 EXP. DATE: 04/2013
 2. RPR REACTIVE (1:4) CONTROL LOT #: 1582 EXP. DATE: 08/2013
 3. RPR NON-REACTIVE CONTROL LOT #: 4371 EXP. DATE: 08/2013


TOLERANCE LEVELS OF CONTROLS

<u>CONTROL</u>	<u>LOT #</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>
REACTIVE	1582	08/2013	REACTIVE
NON-REACTIVE	4371	08/2013	NON-REACTIVE

TEST VALUES/RESULTS

<u>CONTROL</u>	<u>RESULT 1</u>	<u>RESULT 2</u>	<u># PASS</u>
REACTIVE	REACTIVE	REACTIVE	2/2
NON-REACTIVE	NON-REACTIVE	NON-REACTIVE	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	<u>05/04/2011</u>
-----------------	--	-------	-------------------