



Certificate of Analysis

Product Description: RPR (500T)
Code No: SYRPR500
Batch No: RK-1611-6
Manufacture Date: 2016/11
Expiry: 2018/07

CDC Reactive Human Plasma

	1/1	1/2	1/4	1/8	1/16	1/32	1/64
Expected Titre	SR	WR	TR	N	N	N	N
Achieved Titre	SR	WR	TR	N	N	N	N

Specificity

Tested: 50 x Donor Plasma
Results: 50 x Negative

Biohazard

Positive and Negative Controls are produced from pooled human plasma. All human blood components have been tested negative for HBsAg, HCV and HIV 1 & 2 antibodies by methods approved by the FDA. However as no known test methods offer total assurance, materials should be handled as though capable of transmitting disease.

Product Release:

The above batch has passed the quality control specifications and is released for use.

Signed: *Cherith Morris* (QC Manager) **Date:** 2016/11