



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URS-10 CLICHEK**

MANUFACTURE DATE: 02/15/2010

CATALOG#: URS-10-100CK

LOT #: 23784

EXP. DATE: 05/2012

COMPONENTS:	LEUKOCYTES	LOT #:	23223	EXP. DATE:	06/2012
	NITRITE	LOT #:	22863	EXP. DATE:	06/2012
	UROBILINOGEN	LOT #:	22984	EXP. DATE:	06/2012
	PROTEIN	LOT #:	23445	EXP. DATE:	07/2012
	pH	LOT #:	23781	EXP. DATE:	07/2012
	BLOOD	LOT #:	23222	EXP. DATE:	06/2012
	SPECIFIC GRAVITY	LOT #:	22908	EXP. DATE:	06/2012
	KETONE	LOT #:	22518	EXP. DATE:	05/2012
	BILIRUBIN	LOT #:	23257	EXP. DATE:	06/2012
	GLUCOSE	LOT #:	23739	EXP. DATE:	07/2012

TOLERANCE LEVELS OF CONTROLS

Controls	Leuko-cytes	Nitrite	Urobili-nogen	Protein	pH	Blood	Specific Gravity	Ketone	Bilirubin
KOVA I LOT#: 134042 EXP: 08/2010	2+ to 3+	POS (+)	4.0 to $\geq$ 8.0 EU/dl	100 to $\geq$ 300 mg/dl	7.0 to 8.0	2+ to 3+	1.010 to $\geq$ 1.030	40 to $\geq$ 80 mg/dl	2+ to 3+
KOVA II LOT#: 133585 EXP: 05/2010	Trace to 2+	POS (+)	0.2 to 1.0 EU/dl	NEG (-) to 30 mg/dl	6.5 to 8.5	Trace to 2+	$\leq$ 1.005 to 1.015	Trace to $\geq$ 40 mg/dl	1+ to 2+
KOVA III LOT#: 134047 EXP: 05/2011	NEG (-)	NEG (-)	0.2 to 1.0 EU/dl	NEG (-)	5.0 to 6.5	NEG (-)	$\leq$ 1.005 to 1.015	NEG (-)	NEG (-)

Glucose:

Control Level:	Level 1	Level 2	Level 3	Level 4	Level 5
Target (mg/dL):	0	50	100	300	$\geq$ 1000
Range (mg/dL):	0	25 - 75	50 - 150	150 - 400	500 - $\geq$ 1000


TEST VALUES/RESULTS

Controls	Leukocytes	Nitrite	Urobili-nogen	Protein	pH	Blood	Specific Gravity	Ketone	Bilirubin
KOVA I	2+ 2+	+ +	8.0 8.0	300 300	7.5 7.5	3+ 3+	1.015 1.015	80 80	3+ 3+
KOVA II	1+ 1+	+ +	0.2 0.2	Trace Trace	7.5 7.5	1+ 1+	1.010 1.010	40 40	2+ 2+
KOVA III	- -	- -	0.2 0.2	- -	6.5 6.5	- -	1.010 1.010	- -	- -

Glucose:

Test Results:	0	50	100	300	1000
---------------	---	----	-----	-----	------

COMMENTS: VISUAL CHECK OK.

APPROVE D BY:	 Quality Control	DATE :	<u>02/15/2010</u>
------------------	--	-----------	-------------------