



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URS-10 100 Clicheck**

MANUFACTURE DATE: 11/18/2010

Finished Kit

CATALOG#: URS-10-100CK

LOT #: 29192

EXP. DATE: 07/2012

COMPONENTS:	LEUKOCYTES	LOT #:	26151	EXP. DATE:	10/2012
	NITRITE	LOT #:	26153	EXP. DATE:	10/2012
	UROBILINOGEN	LOT #:	26181	EXP. DATE:	10/2012
	PROTEIN	LOT #:	26049	EXP. DATE:	10/2012
	pH	LOT #:	26237	EXP. DATE:	10/2012
	BLOOD	LOT #:	26003	EXP. DATE:	10/2012
	SPECIFIC GRAVITY	LOT #:	26221	EXP. DATE:	10/2012
	KETONE	LOT #:	24902	EXP. DATE:	07/2012
	BILIRUBIN	LOT #:	26159	EXP. DATE:	10/2012
	GLUCOSE	LOT #:	24874	EXP. DATE:	07/2012

TOLERANCE LEVELS OF CONTROLS

Controls	<u>Leuko-</u> <u>cytes</u>	<u>Nitrite</u>	<u>Urobili-</u> <u>nogen</u>	<u>Protein</u>	<u>pH</u>	<u>Blood</u>	<u>Specific</u> <u>Gravity</u>	<u>Ketone</u>	<u>Bilirubin</u>
KOVA I LOT#: 139036-4 EXP: 01/2012	2+ to 3+	POS (+)	4.0 to $\geq$ 8.0 EU/dl	100 to $\geq$ 300 mg/dl	7.0 to 8.0	2+ to 3+	1.010 to $\geq$ 1.030	40 to $\geq$ 80 mg/dl	2+ to 3+
KOVA II LOT#: 138183-2 EXP: 07/2011	Trace to 2+	POS (+)	0.2 to 1.0 EU/dl	NEG (-) to 30 mg/dl	6.5 to 8.5	Trace to 2+	$\leq$ 1.005 to 1.015	Trace to $\geq$ 40 mg/dl	1+ to 2+
KOVA III LOT#: 139316-9 EXP: 08/2012	NEG (-)	NEG (-)	0.2 to 1.0 EU/dl	NEG (-)	5.0 to 6.5	NEG (-)	$\leq$ 1.005 to 1.015	NEG (-)	NEG (-)

Glucose:

Control Level:	Level 1	Level 2	Level 3	Level 4	Level 5
Target (mg/dL):	0	50	100	300	$\geq$ 1000
Range (mg/dL):	0	25 - 75	50 - 150	150 - 400	500 - $\geq$ 1000

TEST VALUES/RESULTS

<u>Controls</u>	<u>Leukocytes</u>	<u>Nitrite</u>	<u>Urobili-</u> <u>nogen</u>	<u>Protein</u>	<u>pH</u>	<u>Blood</u>	<u>Specific</u> <u>Gravity</u>	<u>Ketone</u>	<u>Bilirubin</u>
KOVA I	2+ 2+	+ +	8.0 8.0	300 300	7.5 7.5	3+ 3+	1.020 1.020	2+ 2+	3+ 3+
KOVA II	2+ 2+	+ +	0.2 0.2	Trace Trace	7.0 8.5	1+ 1+	1.015 1.015	1+ 1+	1+ 1+
KOVA III	- -	- -	0.2 0.2	- -	6.0 6.5	- -	1.005 1.005	- -	- -

Glucose:

Test Results:	0	50	100	300	1000
---------------	---	----	-----	-----	------

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	<u>11/18/2010</u>
-----------------	--	-------	-------------------