

TECO DIAGNOSTICS

1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URS-1K 100 Teco**
Finished Kit

MANUFACTURE DATE: 03/22/2011

CATALOG #: URS-1K-100

LOT #: 31535

EXP. DATE: 07/2013

COMPONENT: URS Paper - Ketone

LOT#: 30881

EXP. DATE: 07/2013

TOLERANCE LEVELS OF CONTROLS


CONTROLS	LOT #	EXP. DATE	RANGE
KOVA I	139036-4	01/2012	40 - $\geq$ 160 mg/dL
KOVA II	139317-7	02/2012	Trace - $\geq$ 40 mg/dL
KOVA III	140400-1	03/2013	Negative (-)

TEST VALUES/RESULTS

CONTROLS	1	2	# PASS
KOVA I	80	80	2/2
KOVA II	15	15	2/2
KOVA III	-	-	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:


Quality Control

DATE: 03/22/2011