



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **URS-3 100 Teco Finished Kit** MANUFACTURE DATE: 04/22/2011
CATALOG #: URS-3-100
LOT #: 32173 EXP. DATE: 07/2013

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT #:</u>	<u>EXP. DATE:</u>
URS Paper - Protein	URSP-PRO	31765	08/2013
URS Paper - Glucose	URSP-GLU	31764	08/2013
URS Paper - pH	URSP-PH	31407	07/2013

TOLERANCE LEVELS OF CONTROLS

CONTROLS	LOT #	EXP. DATE	PROTEIN	GLUCOSE	pH
KOVA I	139036-4	01/2012	100 to $\geq$ 2000 mg/dl	500 to $\geq$ 2000 mg/dl	7.0 – 8.5
KOVA II	139317-7	02/2012	Negative (-) to 30 mg/dl	Trace to 500 mg/dl	6.5 – 8.5
KOVA III	140400-1	03/2013	Negative (-)	Negative (-)	5.0 – 6.5

TEST VALUES/RESULTS

CONTROLS	PROTEIN		GLUCOSE		pH		# PASS
	<u>Result 1</u>	<u>Result 2</u>	<u>Result 1</u>	<u>Result 2</u>	<u>Result 1</u>	<u>Result 2</u>	
LEVEL I	300	300	1000	1000	7.5	7.5	6 / 6
LEVEL II	Trace	Trace	100	250	7.0	7.0	6 / 6
LEVEL III	-	-	-	-	6.0	6.0	6 / 6

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:


Quality Control

DATE: 04/22/2011